Vario

- EN** SERVICE MANUAL
- DE** SERVICEANLEITUNG
- FR** INSTRUCTIONS DE SERVICE

Vario 8 (AC/DC)

	Vario 8	Vario 8	Vario 8 c/i	Vario 18	Vario 18	Vario 18 c/i
		AC/DC	AC/DC		AC/DC	AC/DC
 vacuum	low vacuum - 90 mbar - 9 kPa - 68 mmHg	low vacuum - 90 mbar - 9 kPa - 68 mmHg	low vacuum - 90 mbar - 9 kPa - 68 mmHg	high vacuum - 75 kPa - 563 mmHg	high vacuum - 75 kPa - 563 mmHg	medium vacuum - 55 kPa - 413 mmHg
 flow	low flow 8 l/min.	low flow 8 l/min.	low flow 8 l/min.	low flow 18 l/min.	low flow 18 l/min.	low flow 18 l/min.

AC-Version

Vario 8 / Vario 18		230 – 240 V	50 Hz	90 W		T 0,8 A	Ø 5 x 20 mm
Vario 8 / Vario 18		230 – 240 V	60 Hz	90 W		T 0,8 A	Ø 5 x 20 mm
Vario 8 / Vario 18		120 V	60 Hz	70 W		T 1,25 A	Ø 6,3 x 32 mm

AC/DC-Version

Vario 8 / Vario 18		100 – 240 V	50/60 Hz	80 VA		T 1,0 A	Ø 5 x 20 mm
Vario 8 c/i / Vario 18 c/i							

 +5 °C to +40 °C	Operation Betrieb Fonctionnement	 70 °C to 106 °C	Operation Betrieb Fonctionnement
 -20 °C to +50 °C	Transport/Storage Transport/Lagerung Transport/Stockage	 20 °C to 95 °C	Transport/Storage Transport/Lagerung Transport/Stockage
 380 x 170 x 285 mm		 Vario 8: 49 dB(A) / Vario 18: 52 dB(A) Max. noise level Max. Lautstärke Niveau de bruit maximum	
 AC 3.5 kg AC/DC 4.2 kg		 CLASS	 IPX1
		 ISO 9001 ISO 13485 CE (93/42/EEC), IIa	 CE 0123

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Vacuum gauge
Manometer
Manomètre

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I Important notes

In the application area of medical suction devices, the vacuum is normally quoted as the difference (in absolute figures) between the absolute pressure and the atmospheric pressure or as negative values in kilopascals (kPa). The indication, e.g. of -10 kPa, in this document therefore always refers to a pressure level in kPa below the ambient atmospheric pressure.

Rechargeable battery: charge the battery before use. -----

II Use in accordance with the instructions

The device is highly suitable for the aspiration of surgical fluids, secretions, tissue (including bone) and infectious particles. In the emergency room, the doctor's practice, the ward or in the home care area.

	Vario 8	Vario 8 AC/DC	Vario 8 c/i*	Vario 18	Vario 18 AC/DC	Vario 18 c/i*
General suction purposes				●	●	●
Thoracic drainage	●	●	●			
Wound drainage				●	●	●
ENT indications				●	●	●
Gastric suction						●

*c = constant suction / i = intermittent suction

An extensive range of accessories makes the Vario the ideal device for a broad range of medical applications. Contact us – we will be pleased to advise you.

III Guarantee

2 years guarantee on the device from delivery date – 6 months on rechargeable battery and Quatroflex™ suction unit. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage as well as use by unauthorized persons.

IV Service work

- Service work is to be performed only by Medela authorized personnel
- Medela recommendation: perform a routine check once every year -----

V Cleaning

■ Before cleaning the device, remove the mains plug from the socket!

- Comply with in-house directives on hygiene
- Clean the device with a damp cloth and with non-aggressive cleaning material
- Suction jars, lids, suction bags and other accessories are to be handled in accordance with the instruction-sheets provided
- The Medela cleaning recommendations must also be respected

077.0225

VI Functional test

- a. Plug device in and switch on

With c/i versions,
select continuous vacuum

- b. Set maximum vacuum

Press knob to turn

- c. Close vacuum port

- d. Compare vacuum

			low vacuum	medium vacuum	high vacuum
			M I N I M	U M	V A L U E S
+ 2000 m			795 hPa	- 7,0 kPa	- 43 kPa
+ 1000 m			899 hPa	- 7,9 kPa	- 48 kPa
+ 500 m			955 hPa	- 8,5 kPa	- 51 kPa
0 m			1013 hPa	- 9,0 kPa	- 55 kPa
				- 75 kPa	

ok? ▶

Minimum vacuum
values not reached

VII Test filter

Filter open!

1. Set maximum vacuum
2. Switch on the Vario
3. Compare vacuum

low vacuum	medium vacuum	high vacuum
L I	M I	T S
-3 kPa	-10 kPa	-10 kPa

ok? ▶

Vacuum values exceed limits,
replace filter

077.0571

Overflow protection / bacteria filter

VIII Routine check

(time required approx. 30 minutes)

1. Functional test

- Is motor running?
- Is vacuum sufficient?

ok?

▶ not OK

Motor not running
Vacuum insufficient

2. Open device

- Disconnect the pump from the mains supply
- Release four screws
- Remove rear panel

3. Check QuatroFlex™

- Film hinge torn
- open valve cover

ok?

▶ not OK

Replace QuatroFlex™

077.1700

4. Check sound absorber

- Dirt, discolouration
- Brittle material

ok?

▶ not OK

Replace sound absorber

077.1702

5. Check tubings

- Check seal of all connections
- Dirt, discolouration
- Brittle material

ok?

▶ not OK

Replace sound absorber
Replace QuatroFlex™

077.1702

077.1700

6. Check power cord

- Damage, kinks
- Check for tight seal

ok?

▶ not OK

Replace power cord

7. Close device

- Clean interior of device
- Attach rear panel
- Tighten four screws
- Plug in the mains cord
- Only for AC/DC versions:
Testing of rechargeable battery

8. Functional test

IX Battery management

(only for AC/DC versions)

1. Performance of rechargeable battery

	low vacuum	medium vacuum	high vacuum
Duration of pumping	~ 120 min.	~ 30 min.	~ 30 min.
Charging time (charged up to 100%)	~ 5 hrs.	~ 5 hrs.	~ 5 hrs.

2. Testing of rechargeable battery

	low vacuum	medium vacuum	high vacuum
1. Disconnect the mains cord			
2. Switch on the Vario			
3. Plug the suction port	-90 mbar - 9 kPa	-55 kPa	-55 kPa
4. Set required vacuum	-68 mmHg	-413 mmHg	-413 mmHg
5. Measure the running time			
6. Compare the running time	> 120 min.	> 30 min.	> 30 min.

ok?

OK

if above values are reached

Do not forget to charge battery before use

not OK

Conditioning the rechargeable battery

- A Completely discharge the rechargeable battery
- Disconnect the mains cord
 - Switch on the Vario
 - Let the Vario run to a standstill
- B Charge battery
- Plug in the mains cord
 - the rechargeable battery is being charged
 - Wait until the charging indicator goes out
- C Test rechargeable battery

ok?

OK

If above values are reached

Do not forget to charge battery before use

not OK, change rechargeable battery

A1

Storage of rechargeable batteries

In order to keep the self-discharge at a minimum, store the Vario and the replacement rechargeable batteries at temperatures below 25°C (77°F).
Repeat Steps A and B every 60-90 days.

X Insufficient vacuum

1. Open device

2. Check tubings

- Tubing connections
- Dirt, discolouration
- Brittle material

ok?

not OK

Attach tubings

Replace sound absorbers
Replacing QuatroFlex™

077.1702

077.1700

3. Check vacuum port

- Brittle material
- Vacuum port broken

ok?

not OK

Replace vacuum port

077.1711

4. Check QuatroFlex™

- Film hinge torn
- Open valve cover

ok?

not OK

- Disconnect 3 tubings
- Turn QuatroFlex™ 45° counter-clockwise and remove
- General inspection and cleaning

d. New QuatroFlex™

077.1700

- Position on drive shaft
- Press down against holder
- Turn 45° clockwise as far as limit/re-attach all the tubings

5. Check vacuum gauge/regulator

- Mechanical damage
- Oxidation/rust

ok?

not OK

Replace vacuum gauge
Replace regulator

6. Close device

7. Functional test

XI Motor not running

1. Check fuse

- Fuse fitted
- Fuse defective

ok?

not OK →

Fit correct fuse as per technical specifications. See technical data (envelope)

2. Open device

VIII,2

3. Dismantle drive unit

- Dismantle QuatroFlex™
- Dismantle sound absorber
- AC/DC version: disconnect battery
- Dismantle drive unit with 4 screws

4. Check mains cord

- Mains voltage correct
- Correctly connected
- Broken cable
- Electrical connections

ok?

not OK →

Tighten electrical connections

Replace mains cord

A1

5. Check motor/drive belt

- Motor correctly connected
- Motor defective
- Drive belt positioned correctly
- Drive belt torn

ok?

not OK →

Plug in motor
Position drive belt

Replace motor

Replace drive belt

A1

077.1703

6. On/Off switch

- Correctly connected
- Defective

ok?

not OK →

Mount switch

Replace switch

A1

Continued on next page

Motor not running (continued)

AC version
only

- 7. Check capacitor**
- Correctly connected?

ok?

not OK

Connect cord

Replace capacitor

AC/DC version
only

- 8. Check electronic module**
- Correctly connected
- Defective
- Check the jumper setting

ok?

not OK

Connect cable properly
Set jumper correctly

Replace defective
electronic module

077.1765

9. Mount drive unit

- secure drive unit with 4 screws
- AC/DC version: mount battery
- mount QuatroFlex™
- Mount sound absorber

10. Close device

11. Functional test


XII Jumper settings

- Check the jumper setting. Also see the jumper settings (foldout)
- Check correct position of the On/Off switch

c/i-versions

On	Off
8 sec.	4 sec.
* 16 sec.	8 sec.
32 sec.	16 sec.

Vario 8	Vario 18
* with	without

* Factory setting

XIII Overflow

The following service work must be performed immediately to prevent a possible infection. Comply with Medela hygiene instructions

1. Open device

2. Remove and replace following components:

- QuatroFlex™
- Regulator
- Vacuum gauge
- Sound absorber
- Vacuum port

3. General cleaning

Clean device on inside and outside

4. Close device

5. Functional test

XIV Disposal

Before disposal, the device is to be rendered unusable and it must be disposed of in accordance with the local regulations. Must not be disposed of together with household refuse. Hand the discarded device to a local collection and recycling centre for old electrical and electronic appliances. Incorrect disposal can have harmful effects on the environment and public health.

XV Safety-related checks

For the safety-related checks, it is assumed that the device is maintained and repaired throughout its service life strictly and exclusively in compliance with these service instructions by authorized service centres

For the Vario as a device in protection class II (EN IEC 60601-1, A1, A2:1995), the safety-related checks are confined to visual inspection of the housing and mains power cord for damage. This check must be performed prior to each use.

Devices of protection class II do not have a protective earth conductor; there is therefore no need to check the earth leakage current!

The Vario enclosure is made entirely of insulating material. Tests of the **enclosure leakage current** using common measuring instruments will therefore not reveal measurable values.

The applied parts of the Vario (type B or CF) are connected to the Vario via non-conductive vacuum hoses, collection jars and intermediate tubings.

Even when suctioning a conductive fluid until the overflow protection device activates (9 g/l NaCl; EN ISO 10079-1;1999), measurements of the **patient leakage current** using common measuring instruments will not reveal measurable values.

The Vario pumps do not have patient circuits, signal input and output parts and functional earth connections.

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Notes • Notizen • Notes

A3 Occupational Safety and Health

Guideline for Service Workshops of Medela and Representatives

The biggest problem with any type of disease-causing organism is that it cannot be seen with the naked eye. Some people believe that, if we cannot see the disease-causing organism, it will not cause us any problems. THIS IS FALSE!

There are numerous Medela medical devices in circulation: surgical aspirators, home care pumps, breast pumps, etc. The possible exposition of Medela employees to contamination is therefore imminent.

We strongly recommend workshop personnel and other staff working on potentially contaminated equipment to be vaccinated against hepatitis B.
Carefully consider further protection provided by vaccination and other precautionary measures which may be advisable in your region.

When devices are returned to Service Shops for repairs, the device/components are to be considered hazardous and potentially infectious!

Such equipment shall be separated, marked accordingly and be stored in an isolated location until disinfection.

Disinfection Procedure:

- Always wear suitable work clothing, protective gloves and goggles
- Unpack equipment and place it on a suitable bench (sink)
- Open device and dispose of any soiled, contaminated parts like internal tubing, valve plates, etc. according to local regulations
- Disassemble piston/cylinder units which are not worn out.

All possibly contaminated parts:

1. Spray the surfaces with disinfectant
(Caltech Citrex®, ALDA Pharmaceuticals Viralex™, Schülke & Mayr Mikrocid® AF Liquid or similar Hospital Spray Disinfectants) and leave for 3 minutes
 2. Clean the surfaces with a good cleaning agent, methodically rinse with water
 3. Again spray the exterior with disinfectant and wipe parts dry with household tissue or allow to air dry
- Dispose of protective gloves after the disinfection procedure
 - Wash your hands

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Medela UK Ltd.
Medela, Inc.
Medela AG, Medical Technology